

合成吸収性縫合糸

LACLON

ラクロン

販売名：P(LA/CL) 縫合糸
高度管理医療機器 クラスⅣ
承認番号：22800BZX00095000
色：紫、クリア
号数：EP4～EP0.5
(USP1～USP7-0)





















丸針 TAPER POINT



針サイズ	入数	色	糸長	糸の太さ E.P. (U.S.P.)																
				0.5 (7-0)	0.7 (6-0)	1 (5-0)	1.5 (4-0)	2 (3-0)	3 (2-0)	3.5 (1-0)										
D 弱弯 3/8 Circle																				
D-8mm	両端	1本×10袋	紫	60cm	LD507															
D-9mm	両端	1本×30袋	紫	60cm	LD527	LD516														
				75cm	LD517	LD506														
D-11mm	両端	1本×10袋	紫	30cm				LD504												
				75cm		LD526														
D-13mm	両端	1本×10袋	紫	60cm			LD515													
				1本×30袋	60cm		LD546													
					75cm		LD536	LD505												
D-17mm		1本×30袋	紫	70cm				LD014												
C 強弯 1/2 Circle																				
C-10mm	両端	1本×20袋	紫	70cm		LC506														
C-11mm		1本×30袋	紫	70cm			LC025													
C-13mm		1本×30袋	紫	70cm		LC026		LC044												
				8本×10袋 トリプル針	45cm			LCT005	LCT024											
C-13mm	両端	1本×30袋	紫	70cm			LC515													
				75cm		LC516														
C-17mm		1本×10袋	紫	60cm			LC045													
				1本×30袋	70cm			LC035	LC014	LC003										
					8本×10袋 トリプル針	45cm			LCT015	LCT014	LCT003									
				60cm				LCT025	LCT004	LCT023										
				70cm						LCT083										
C-17mm	両端	1本×30袋	紫	75cm				LC525												
				90cm			LC535	LC514												

※上記規格品以外は、担当者へお問い合わせ下さい。





■ 丸針 TAPER POINT

針サイズ	入数	色	糸長	糸の太さ E.P. (U.S.P.)						
				0.7 (6-0)	1 (5-0)	1.5 (4-0)	2 (3-0)	3 (2-0)	3.5 (1-0)	4 (1号)
C 強弯 1/2 Circle										
C-22mm 	1本×30袋	紫	70cm			LC024	LC013			
	8本×10袋 	紫	45cm			LCT044	LCT013			
		紫	60cm			LCT054	LCT063			
C-22mm  両端	1本×30袋	紫	90cm			LC534	LC513			
C-26mm 	1本×30袋	紫	70cm			LC034	LC023	LC002		
	8本×10袋 	紫	45cm			LCT064	LCT033	LCT002		
		紫	70cm				LCT073			
C-26mm  両端	1本×30袋	紫	70cm				LC543			
C-26mm  両端	1本×30袋	紫	90cm				LC533			
C-30mm 	1本×30袋	紫	70cm			LC033	LC012			
C-36mm 	1本×30袋	紫	70cm			LC043	LC022	LC001	LC000	
	8本×10袋 	紫	45cm				LCT053	LCT012	LCT001	LCT000
C-36mm 	1本×10袋	紫	150cm					LC022L		
C-40mm 	1本×30袋	紫	70cm					LC011	LC010	
	8本×10袋 	紫	45cm					LCT011	LCT010	
C-40mm 	1本×20袋	紫	150cm					LC011L		
C-48mm 	1本×20袋	紫	70cm						LC020	
	8本×10袋 	紫	45cm					LCT021	LCT020	
C-48mm 	1本×20袋	紫	90cm					LC021L	LC020L	
		紫	150cm					LC031L		
C-65mm 	1本×10袋	紫	240cm						LC030L	

※上記規格品以外は、担当者へお問い合わせ下さい。




ブランテーパー BLUN TAPER



針サイズ	入数	色	糸長	糸の太さ E.P. (U.S.P.)							
				0.7 (6-0)	1 (5-0)	1.5 (4-0)	2 (3-0)	3 (2-0)	3.5 (1-0)	4 (1号)	
GC 強弯 1/2 Circle											
GC-26mm 	8本×10袋 トローリング針	紫	45cm					LGCT003	LGCT002		
			70cm							LGCT001	
GC-36mm 	1本×30袋	紫	90cm						LGCT022	LGCT041	LGCT000
			8本×10袋 トローリング針	45cm						LGCT012	LGCT001
GC-40mm 	8本×10袋 トローリング針	紫	45cm							LGCT011	LGCT010
GC-48mm 	1本×20袋	紫	90cm							LGCT021	LGCT010
			8本×10袋 トローリング針	45cm							LGCT021
















鈍針 BLUNT POINT



針サイズ	入数	色	糸長	糸の太さ E.P. (U.S.P.)							
				0.7 (6-0)	1 (5-0)	1.5 (4-0)	2 (3-0)	3 (2-0)	3.5 (1-0)	4 (1号)	
G 強弯 1/2 Circle											
G-48mm 	1本×10袋	紫	150cm							LG031L	
			240cm								
G-65mm 	1本×10袋	紫	90cm							LG001	
G-65mm 	1本×10袋	紫	240cm								LG000L

角針 REVERSE CUTTING



針サイズ	入数	色	糸長	糸の太さ E.P. (U.S.P.)							
				0.5 (7-0)	0.7 (6-0)	1 (5-0)	1.5 (4-0)	2 (3-0)	3 (2-0)	3.5 (1-0)	4 (1号)
B 弱弯 3/8 Circle											
B-11mm 	1本×10袋	クリア	45cm		LB006						
	1本×20袋	クリア	45cm			LB005					
B-13mm 	1本×10袋	クリア	45cm			LB025	LB014				
	1本×30袋	紫	45cm			LB015	LB004				
B-15mm 	1本×20袋	クリア	45cm			LB035	LB024				
B-20mm 	1本×10袋	紫	45cm				LB034				
	1本×20袋	クリア	45cm				LB054				
	1本×30袋	クリア	70cm				LB044	LB003			
B-24mm 	1本×30袋	クリア	70cm				LB064				
A 強弯 1/2 Circle											
A-8mm 	1本×10袋	紫	45cm	LA007	LA006						
		クリア	45cm		LA016	LA005	LA004				
A-10mm 	1本×10袋	クリア	45cm		LA026	LA015					
A-13mm 	1本×10袋	紫	45cm			LA025	LA014				
		クリア	45cm			LA035	LA024				
	1本×30袋	クリア	45cm		LA036						
A-15mm 	1本×10袋	クリア	45cm			LA045	LA034	LA003			
A-22mm 	1本×20袋	クリア	70cm					LA013			
	8本×10袋 	紫	45cm				LAT004	LAT003			
A-24mm 	1本×10袋	クリア	70cm					LA023	LA002		
A-26mm 	1本×30袋	紫	70cm					LA033	LA012		
	8本×10袋 	紫	45cm					LAT013	LAT002		
A-36mm 	1本×30袋	紫	70cm							LA001	